

## CAB24440

---

### Product Information

<b>Product SKU:</b>	CAB24440	<b>Gene ID:</b>	23308	<b>Size:</b>	100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Mouse

---

### Additional Information

<b>Observed MW:</b>	65-70kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	33kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Enables identical protein binding activity. Predicted to be involved in T cell receptor signaling pathway and positive regulation of interleukin-4 production. Located in cytoplasmic ribonucleoprotein granule and plasma membrane.

**Recommended Dilution:** WB, 1:500 - 1:1000

**Synonyms:** B7h; B7H2; GL50; B7-H2; B7RP1; CD275; ICOSL; LICOS; B7RP-1; ICOS-L

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 19-256 of human ICOSL (NP\_056074.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.