

CAB24458

Product Information

| | | |
|------------------------------|-----------------------------|--------------------------|
| Product SKU: CAB24458 | Gene ID: - | Size: 20uL, 100uL |
| Clone No: - | Host Species: Rabbit | Reactivity: Human |

Additional Information

| | |
|----------------------------|--------------------------------|
| Observed MW: >10kDa | Conjugate: Unconjugated |
| Calculated MW: - | Isotype: IgG |

Immunogen Information

Background:

Recommended Dilution: WB, 1:100 - 1:500

Synonyms: -

Purification Method: Affinity purification

Immunogen: A synthesized peptide derived from plxIS motif.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.