

CAB24465

Product Information

Product SKU: CAB24465	Gene ID: -	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human

Additional Information

Observed MW: 30kDa	Conjugate: Unconjugated
Calculated MW: -	Isotype: IgG

Immunogen Information

Background: S tag is an epitope tag composed of a 15 residue peptide, KETAAAKFERQHMD_S, derived from the pancreatic ribonuclease A.

Recommended Dilution: WB, 1:1000 - 1:5000

Synonyms: S-Tag; S

Purification Method: Affinity purification

Immunogen: A synthesized peptide derived from S-Tag.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.