

## CAB24468

---

### Product Information

<b>Product SKU:</b>	CAB24468	<b>Gene ID:</b>	79977	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Rat

---

### Additional Information

<b>Observed MW:</b>	71kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	69kDa/71kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** The protein encoded by this gene is a transcription factor that can act as a homodimer or as a heterodimer with either GRHL1 or GRHL3. Defects in this gene are a cause of non-syndromic sensorineural deafness autosomal dominant type 28 (DFNA28).

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** GRHL2; BOM; DFNA28; ECTDS; TFCP2L3; grainyhead-like protein 2 homolog

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 110-210 of human GRHL2.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.