

CAB24470

Product Information

Product SKU:	CAB24470	Gene ID:	84913	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	35kDa	Conjugate:	Unconjugated
Calculated MW:	35kDa	Isotype:	IgG

Immunogen Information

Background:

Recommended Dilution: WB, 1:1000 - 1:5000

Synonyms: ATOH8; HATH6; bHLHa21; protein atonal homolog 8

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ATOH8 (NP_116216.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.