

CAB24471

Product Information

Product SKU: CAB24471	Gene ID: 65108	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 40kDa	Conjugate: Unconjugated
Calculated MW: 19kDa	Isotype: IgG

Immunogen Information

Background: This gene encodes a member of the myristoylated alanine-rich C-kinase substrate (MARCKS) family. Members of this family play a role in cytoskeletal regulation, protein kinase C signaling and calmodulin signaling. The encoded protein affects the formation of adherens junction. Alternative splicing results in multiple transcript variants. Pseudogenes of this gene are located on the long arm of chromosomes 6 and 10.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: MARCKSL1; F52; MACMARCKS; MLP; MLP1; MRP; MARCKS like 1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 95-195 of human MARCKSL1 (NP_075385.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.