

CAB24473

Product Information

| | | | | | |
|---------------------|----------|----------------------|--------|--------------------|-------------|
| Product SKU: | CAB24473 | Gene ID: | 158866 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human |

Additional Information

| | | | |
|-----------------------|-------------------|-------------------|--------------|
| Observed MW: | 36kDa | Conjugate: | Unconjugated |
| Calculated MW: | 17kDa/38kDa/39kDa | Isotype: | IgG |

Immunogen Information

Background: The protein encoded by this gene belongs to the DHHC palmitoyltransferase family. Mutations in this gene are associated with mental retardatio X-linked type 91 (MRX91). Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: ZDHHC15; DHHC15; MRX91; palmitoyltransferase ZDHHC15

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 231-331 of human ZDHHC15 (NP_659406.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.