

CAB24474

Product Information

Product SKU:	CAB24474	Gene ID:	55742	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	45-50kDa	Conjugate:	Unconjugated
Calculated MW:	20kDa/42kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the parvin family of actin-binding proteins. Parvins are associated with focal contacts and contain calponin homology domains that bind to actin filaments. The encoded protein is part of the integrin-linked kinase signaling complex and plays a role in cell adhesion, motility and survival.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: PARVA; CH-ILKBP; MXRA2; alpha-parvin

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PARVA (NP_060692.3) .

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.