

CAB24477

Product Information

Product SKU: CAB24477	Gene ID: 207/208/10000	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 60kDa	Conjugate: Unconjugated
Calculated MW: -	Isotype: IgG

Immunogen Information

Background:

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: -

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-123 of human AKT1/AKT2/AKT3 (NP_005154.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.