

## CAB24489

---

### Product Information

<b>Product SKU:</b>	CAB24489	<b>Gene ID:</b>	17123	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Mouse

---

### Additional Information

<b>Observed MW:</b>	60kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	43kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Predicted to enable integrin binding activity involved in cell-matrix adhesion. Acts upstream of or within keratinocyte differentiation and leukocyte migration. Predicted to be located in membrane. Predicted to be integral component of membrane. Is expressed in eyelid; hair follicle; and hemolymphoid system. Orthologous to human MAdCAM1 (mucosal vascular addressin cell adhesion molecule 1).

**Recommended Dilution:** WB, 1:500 - 1:1000

**Synonyms:** MAdCAM-1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 22-109 of mouse MAdCAM-1 (NP\_038619.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.