

CAB24493

Product Information

| | | | | | |
|---------------------|----------|----------------------|--------|--------------------|-------------|
| Product SKU: | CAB24493 | Gene ID: | 171283 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human |

Additional Information

| | | | |
|-----------------------|-------|-------------------|--------------|
| Observed MW: | 50kDa | Conjugate: | Unconjugated |
| Calculated MW: | 33kDa | Isotype: | IgG |

Immunogen Information

Background: Enables phosphatidylserine binding activity. Acts upstream of or within positive regulation of mast cell activation; response to lipopolysaccharide; and response to wounding. Located in brush border and cell surface. Is expressed in brain; dorsal aorta; liver; lung; and reproductive system. Human ortholog(s) of this gene implicated in atopic dermatitis. Orthologous to human HAVCR1 (hepatitis A virus cellular receptor 1).

Recommended Dilution: WB, 1:2000 - 1:8000

Synonyms: Tim1; KIM-1; TIM-1; Timd1; TIM-1/KIM-1/HAVCR

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 22-237 of mouse TIM-1/KIM-1/HAVCR (NP_599009.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.