

CAB24502

Product Information

Product SKU: CAB24502	Gene ID: 9921	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 95kDa	Conjugate: Unconjugated
Calculated MW: 89kDa/90kDa	Isotype: IgG

Immunogen Information

Background: The protein encoded by this gene contains a ring finger motif, which is known to be involved in protein-protein interactions. The specific function of this protein has not yet been determined. EST data suggests the existence of multiple alternatively spliced transcript variants, however, their full length nature is not known.

Recommended Dilution: WB,1:1000 - 1:5000 IF/ICC,1:50 - 1:200

Synonyms: RNF10; RIE2; RING finger protein 10

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 660-760 of human RNF10 (NP_055683.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.