

CAB24504

Product Information

Product SKU:	CAB24504	Gene ID:	25419	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	27kDa	Conjugate:	Unconjugated
Calculated MW:	25kDa	Isotype:	IgG

Immunogen Information

Background:	Enables cholesterol binding activity; identical protein binding activity; and low-density lipoprotein particle binding activity. Involved in several processes, including cellular response to nitric oxide; complement activation, classical pathway; and negative regulation of superoxide anion generation. Located in extracellular space; filopodium; and growth cone. Used to study rheumatoid arthritis. Biomarker of several diseases, including acute kidney failure; artery disease (multiple); iron deficiency anemia; obesity; and status epilepticus. Human ortholog(s) of this gene implicated in several diseases, including Kawasaki disease; autoimmune disease (multiple); macular degeneration (multiple); middle cerebral artery infarction; and nephritis (multiple). Orthologous to human CRP (C-reactive protein).
Recommended Dilution:	WB,1:1000 - 1:5000
Synonyms:	Aa1249; Ac1262; Ab1-341; Ab2-196; Ac1-114; Ac2-069; Ba2-693; Rat C-Reactive Protein/CRP
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 20-230 of rat C-Reactive Protein/CRP (NP_058792.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.