

CAB24505

Product Information

Product SKU:	CAB24505	Gene ID:	23067	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	208kDa/212kDa	Isotype:	IgG

Immunogen Information

Background: SET1B is a component of a histone methyltransferase complex that produces trimethylated histone H3 at Lys4 (Lee et al., 2007 [PubMed 17355966]).[supplied by OMIM, Mar 2008]

Recommended Dilution: IHC-P,1:50 - 1:200

Synonyms: SETD1B; KMT2G; Set1B; SET domain containing 1B

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 150-250 of human SETD1B (NP_001340274.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.