

CAB24509

Product Information

Product SKU:	CAB24509	Gene ID:	28959	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	25kDa	Conjugate:	Unconjugated
Calculated MW:	25kDa/29kDa	Isotype:	IgG

Immunogen Information

Background: TMEM176B (transmembrane protein 176B) is a protein-coding gene. An important paralog of this gene is TMEM176A.

Recommended Dilution: WB, 1:1000 - 1:5000

Synonyms: TMEM176B; LR8; MS4B2; transmembrane protein 176B

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-65 of human TMEM176B (NP_001094781.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.