

CAB24520

Product Information

Product SKU:	CAB24520	Gene ID:	57604	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	51kDa	Conjugate:	Unconjugated
Calculated MW:	14kDa	Isotype:	IgG

Immunogen Information

Background: Enables tRNA methyltransferase activity. Predicted to be involved in tRNA wobble uridine modification. Predicted to be active in cytoplasm and nucleus.

Recommended Dilution: WB,1:1000 - 1:5000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: KIAA1456; C8orf79; TRM9L

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 150-250 of human KIAA1456 (NP_065895.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.