

CAB24527

Product Information

Product SKU: CAB24527	Gene ID: 12508	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Mouse

Additional Information

Observed MW: 35-45kDa	Conjugate: Unconjugated
Calculated MW: 24kDa	Isotype: IgG

Immunogen Information

Background: Predicted to enable identical protein binding activity. Involved in positive regulation of myoblast fusion. Located in cell-cell junction and plasma membrane. Is expressed in several structures, including alimentary system; brain; genitourinary system; hemolymphoid system gland; and liver and biliary system. Orthologous to human CD53 (CD53 molecule).

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: Ox-44; Tspan25; CD53

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 107-181 of mouse CD53 (NP_031677.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.