

CAB24546

Product Information

Product SKU:	CAB24546	Gene ID:	94234	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	41kDa	Isotype:	IgG

Immunogen Information

Background: FOXQ1 is a member of the FOX gene family, which is characterized by a conserved 110-amino acid DNA-binding motif called the forkhead or winged helix domain. FOX genes are involved in embryonic development, cell cycle regulation, tissue-specific gene expression, cell signaling, and tumorigenesis.

Recommended Dilution: IHC-P, 1:50 - 1:200

Synonyms: FOXQ1; HFH1; forkhead box Q1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOXQ1 (NP_150285.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.