

CAB24548

Product Information

Product SKU:	CAB24548	Gene ID:	231386	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	110kDa	Conjugate:	Unconjugated
Calculated MW:	85kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable N6-methyladenosine-containing RNA binding activity and mRNA binding activity. Involved in posttranscriptional regulation of gene expression. Acts upstream of or within gamete generation; in utero embryonic development; and regulation of mRNA processing. Located in nucleus. Orthologous to human YTHDC1 (YTH domain containing 1).

Recommended Dilution: WB, 1:1000 - 1:5000

Synonyms: mKIAA1966; A730098D12Rik; Mouse Ythdc1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 597-736 of mouse Ythdc1 (NP_808348.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.