

CAB24660

Product Information

Product SKU:	CAB24660	Gene ID:	27161	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	95kDa	Conjugate:	Unconjugated
Calculated MW:	97kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the Argonaute family of proteins which play a role in RNA interference. The encoded protein is highly basic, and contains a PAZ domain and a PIWI domain. It may interact with dicer1 and play a role in short-interfering-RNA-mediated gene silencing. Multiple transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution: WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200

Synonyms: PPD; Q10; CASC7; EIF2C2; LESKRES; LINC00980; Argonaute-2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 60-160 of human Argonaute-2 (NP_036286.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.