

CAB24732

Product Information

Product SKU:	CAB24732	Gene ID:	5118	Size:	20uL, 100uL
Clone No:	ARC63882	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	48kDa	Conjugate:	Unconjugated
Calculated MW:	48kDa	Isotype:	IgG

Immunogen Information

Background: Fibrillar collagen types I-III are synthesized as precursor molecules known as procollagens. These precursors contain amino- and carboxyl-terminal peptide extensions known as N- and C-propeptides, respectively, which are cleaved, upon secretion of procollagen from the cell, to yield the mature triple helical, highly structured fibrils. This gene encodes a glycoprotein which binds and drives the enzymatic cleavage of type I procollagen and heightens C-proteinase activity.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: PCPE; PCPE1; PCPE-1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 304-449 of human TLR4 (NP_002584.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.