

CAB24844

Product Information

Product SKU:	CAB24844	Gene ID:	3690	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	97kDa/130kDa	Conjugate:	Unconjugated
Calculated MW:	86kDa/87kDa	Isotype:	IgG

Immunogen Information

Background: The ITGB3 protein product is the integrin beta chain beta 3. Integrins are integral cell-surface proteins composed of an alpha chain and a beta chain. A given chain may combine with multiple partners resulting in different integrins. Integrin beta 3 is found along with the alpha IIb chain in platelets. Integrins are known to participate in cell adhesion as well as cell-surface mediated signalling.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: ITGB3; BDPLT16; BDPLT2; CD61; GP3A; GPIIIa; GT; integrin beta-3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 21-718 of human ITGB3(NP_000203.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.