

CAB24854

Product Information

Product SKU:	CAB24854	Gene ID:	6512	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	70kDa	Conjugate:	Unconjugated
Calculated MW:	17kDa/60kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable anion transmembrane transporter activity. Involved in neurotransmitter reuptake. Predicted to be located in plasma membrane. Predicted to be active in glutamatergic synapse. Predicted to be integral component of postsynaptic membrane and integral component of presynaptic membrane.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: SLC1A7; AAAT; EAAT5; solute carrier family 1 member 7

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 480-560 of human SLC1A7 (NP_006662.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.