

CAB24856

Product Information

Product SKU:	CAB24856	Gene ID:	10672	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	33kDa/44kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable D5 dopamine receptor binding activity; G-protein beta/gamma-subunit complex binding activity; and GTPase activity. Predicted to be involved in several processes, including Rho protein signal transduction; activation of phospholipase D activity; and multicellular organism aging. Predicted to act upstream of or within several processes, including branching involved in blood vessel morphogenesis; negative regulation of vascular associated smooth muscle cell migration; and negative regulation of vascular associated smooth muscle cell proliferation. Located in cytosol and nucleus.

Recommended Dilution: IHC-P,1:50 - 1:200

Synonyms: GNA13; G13; G protein subunit alpha 13

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-377 of human GNA13 (NP_006563.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.