

CAB24862

Product Information

| | | | | | |
|---------------------|----------|----------------------|--------|--------------------|-------------|
| Product SKU: | CAB24862 | Gene ID: | 54567 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human |

Additional Information

| | | | |
|-----------------------|-------|-------------------|--------------|
| Observed MW: | 90kDa | Conjugate: | Unconjugated |
| Calculated MW: | 74kDa | Isotype: | IgG |

Immunogen Information

Background: This gene is a homolog of the Drosophila delta gene. The delta gene family encodes Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: DLL4; AOS6; hdelta2; delta-like protein 4

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 27-524 of human DLL4 (NP_061947.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.