

## CAB24873

---

### Product Information

<b>Product SKU:</b> CAB24873	<b>Gene ID:</b> 54532	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Human, Mouse

---

### Additional Information

<b>Observed MW:</b> 120kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 120kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:**

**Recommended Dilution:** WB, 1:500 - 1:1000 IF/ICC, 1:50 - 1:200

**Synonyms:** USP53; ubiquitin specific peptidase 53

**Purification Method:** Affinity purification

**Immunogen:** Recombinant Protein corresponding to a sequence within amino acids 495-582 of human USP53(NP\_061923.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.