

CAB2540

Product Information

Product SKU: CAB2540	Gene ID: 3161	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 90kDa	Conjugate: Unconjugated
Calculated MW: 84kDa	Isotype: IgG

Immunogen Information

Background: The protein encoded by this gene is involved in cell motility. It is expressed in breast tissue and together with other proteins, it forms a complex with BRCA1 and BRCA2, thus is potentially associated with higher risk of breast cancer. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: CD168; IHABP; RHAMM; CD168/RHAMM

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 341-440 of human CD168/RHAMM (NP_001136028.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.