

UPK3B Rabbit Polyclonal Antibody



CAB2629

Product Information

Size:

50uL, 100uL, 200uL

Observed MW:

35kDa

Calculated MW:

28kDa/29kDa/33kDa

Applications:

WB IHC

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:25
- 1:50

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

80761

Uniprot

Synonyms:

UPK3B


Immunogen:

Recombinant protein of human UPK3B


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of 293T cells, using UPK3B antibody (CAB2629) at dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry of paraffin-embedded Human brain using UPK3B antibody (CAB2629) at dilution of 1:25 (20x lens).