

LGALS2 Rabbit Polyclonal Antibody



CAB2734

Product Information

Size:

50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

14kDa

Applications:

WB IHC

Reactivity:

Human

Protein Background

The protein encoded by this gene is a soluble beta-galactoside binding lectin. The encoded protein is found as a homodimer and can bind to lymphotoxin-alpha. A single nucleotide polymorphism in an intron of this gene can alter the transcriptional level of the protein, with a resultant increased risk of myocardial infarction.

Immunogen information

Gene ID:

3957

Uniprot

P05162

Synonyms:

LGALS2; HL14

Antibody Information

Recommended dilutions:

WB 1:200 - 1:1000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant protein of human LGALS2

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
