

# AKAP8L Rabbit Polyclonal Antibody



CAB3050

---

## Product Information

### Size:

50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

65kDa/71kDa

### Applications:

WB IHC

### Reactivity:

Human, Mouse

## Antibody Information

### Recommended dilutions:

WB 1:200 - 1:1000 IHC 1:25  
- 1:100

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

26993

#### Uniprot

Q9ULX6

#### Synonyms:

AKAP8L; HA95; HAP95; NAKAP; NAKAP95

#### Immunogen:

A synthetic Peptide of human AKAP8L

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---