

LPAR4 Rabbit Polyclonal Antibody



CAB3060

Product Information

Size:

50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

41kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Protein Background

This gene encodes a member of the lysophosphatidic acid receptor family. It may also be related to the P2Y receptors, a family of receptors that bind purine and pyrimidine nucleotides and are coupled to G proteins. The encoded protein may play a role in monocytic differentiation.

Immunogen information

Gene ID:

2846

Uniprot

Q99677

Synonyms:

LPAR4; GPR23; LPA4; P2RY9; P2Y5-LIKE; P2Y9

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:25
- 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

A synthetic peptide of human LPAR4

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
