

TET3 Rabbit Polyclonal Antibody



CAB3141

Product Information

Size:

200uL

Observed MW:

Refer to figures

Calculated MW:

77kDa/166kDa/179kDa

Applications:

WB IP

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:200 - 1:500 IP 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

200424

Uniprot

O43151

Synonyms:

TET3; hCG_40738; TET3

Immunogen:

A synthetic peptide of human TET3

Storage:

Store at 4°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

Product Images
