

CAB3214

Product Information

Product SKU:	CAB3214	Gene ID:	8192	Size:	20uL, 100uL
Clone No:	ARC1929	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	26kDa	Conjugate:	Unconjugated
Calculated MW:	30kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene belongs to the peptidase family S14 and hydrolyzes proteins into small peptides in the presence of ATP and magnesium. The protein is transported into mitochondrial matrix and is associated with the inner mitochondrial membrane.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: DFNB81; PRLTS3; CLPP

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 178-277 of human CLPP (Q16740).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.