

CAB3312

Product Information

Product SKU: CAB3312	Gene ID: 4718	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 14kDa	Conjugate: Unconjugated
Calculated MW: 14kDa	Isotype: IgG

Immunogen Information

Background: Involved in mitochondrial respiratory chain complex I assembly. Located in mitochondrion. Part of mitochondrial respiratory chain complex I. Implicated in nuclear type mitochondrial complex I deficiency.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: HLC-1; B14.5b; MC1DN36; NADHDH2; CI-B14.5b; NDUFC2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NDUFC2 (NP_004540.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.