

CAB3395

Product Information

Product SKU:	CAB3395	Gene ID:	4221	Size:	20uL, 100uL
Clone No:	ARC1968	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	72kDa	Conjugate:	Unconjugated
Calculated MW:	68kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes menin, a tumor suppressor associated with a syndrome known as multiple endocrine neoplasia type 1. Menin is a scaffold protein that functions in histone modification and epigenetic gene regulation. It is thought to regulate several pathways and processes by altering chromatin structure through the modification of histones.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: MEAI; SCG2; Menin

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 516-615 of human Menin (O00255).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.