

CAB3517

Product Information

Product SKU:	CAB3517	Gene ID:	1114	Size:	20uL, 100uL
Clone No:	ARC2023	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	110kDa	Conjugate:	Unconjugated
Calculated MW:	78kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a tyrosine-sulfated secretory protein abundant in peptidergic endocrine cells and neurons. This protein may serve as a precursor for regulatory peptides.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: SCG1; Chromogranin B

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 550-650 of human Chromogranin B (P05060).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.