

## CAB3530

---

### Product Information

<b>Product SKU:</b> CAB3530	<b>Gene ID:</b> 9833	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> ARC2032	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b> 75kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 75kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** Enables calcium ion binding activity; non-membrane spanning protein tyrosine kinase activity; and protein serine/threonine kinase activity. Involved in apoptotic process; cell population proliferation; and protein autophosphorylation. Located in cell cortex and plasma membrane.

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** HPK38; MELK

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MELK (NP\_055606.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.