

CAB3548

Product Information

Product SKU:	CAB3548	Gene ID:	7410	Size:	20uL, 100uL
Clone No:	ARC2041	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	100kDa	Conjugate:	Unconjugated
Calculated MW:	101kDa	Isotype:	IgG

Immunogen Information

Background: VAV2 is the second member of the VAV guanine nucleotide exchange factor family of oncogenes. Unlike VAV1, which is expressed exclusively in hematopoietic cells, VAV2 transcripts were found in most tissues. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: VAV-2; VAV2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human VAV2 (P52735).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.