

CAB3563

---

## Product Information

<b>Product SKU:</b>	CAB3563	<b>Gene ID:</b>	22809	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC2043	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	35kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	31kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** Enables several functions, including DNA-binding transcription activator activity, RNA polymerase II-specific; RNA polymerase II transcription regulatory region sequence-specific DNA binding activity; and tubulin binding activity. Involved in several processes, including fat cell differentiation; regulation of cell cycle process; and regulation of transcription, DNA-templated. Located in centrosome; cytosol; and nucleoplasm.

**Recommended Dilution:** WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

**Synonyms:** ATEFX; HMFN0395; ATF5

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 183-282 of human ATF5 (Q9Y2D1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.