

CAB3566

Product Information

Product SKU:	CAB3566	Gene ID:	4212	Size:	20uL, 100uL
Clone No:	ARC2049	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	50kDa	Conjugate:	Unconjugated
Calculated MW:	52kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a homeobox protein belonging to the TALE ('three amino acid loop extension') family of homeodomain-containing proteins. TALE homeobox proteins are highly conserved transcription regulators, and several members have been shown to be essential contributors to developmental programs. Multiple transcript variants encoding distinct isoforms have been described for this gene.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: MRG1; CPCMR; HsT18361; MEIS2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MEIS2 (O14770).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.