

CAB3579

Product Information

Product SKU:	CAB3579	Gene ID:	4116	Size:	20uL, 100uL
Clone No:	ARC2048	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	17kDa	Conjugate:	Unconjugated
Calculated MW:	17kDa	Isotype:	IgG

Immunogen Information

Background: Drosophila that have mutations in their mago nashi (grandchildless) gene produce progeny with defects in germlasm assembly and germline development. This gene encodes the mammalian mago nashi homolog. In mammals, mRNA expression is not limited to the germ plasm, but is expressed ubiquitously in adult tissues and can be induced by serum stimulation of quiescent fibroblasts.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: MAGOH1; MAGOHA; MAGOH

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 67-146 of human MAGOH (P61326).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.