

CAB3698

Product Information

Product SKU:	CAB3698	Gene ID:	5361	Size:	20uL, 100uL
Clone No:	ARC2077	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	230kDa/250kDa	Conjugate:	Unconjugated
Calculated MW:	211kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable semaphorin receptor activity. Predicted to be involved in several processes, including generation of neurons; regulation of GTPase activity; and regulation of cell shape. Predicted to act upstream of or within dichotomous subdivision of terminal units involved in salivary gland branching; neuron projection morphogenesis; and regulation of smooth muscle cell migration. Located in cytosol and nucleoplasm.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: NOV; NOVp; PLXN1; DWOPNED; PLEXIN-A1; Plexin A1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human Plexin A1 (Q9UIW2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.