

CAB3983

Product Information

Product SKU: CAB3983	Gene ID: 4722	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse

Additional Information

Observed MW: 26kDa	Conjugate: Unconjugated
Calculated MW: 30kDa	Isotype: IgG

Immunogen Information

Background: This gene encodes one of the iron-sulfur protein (IP) components of mitochondrial NADH:ubiquinone oxidoreductase (complex I). Mutations in this gene are associated with Leigh syndrome resulting from mitochondrial complex I deficiency.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: CI-30; MC1DN8; NDUFS3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 37-264 of human NDUFS3 (NP_004542.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.