

CAB4044

Product Information

Product SKU:	CAB4044	Gene ID:	113130	Size:	20uL, 100uL
Clone No:	ARC2096	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	35kDa	Conjugate:	Unconjugated
Calculated MW:	28kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable chromatin binding activity. Involved in double-strand break repair; mitotic sister chromatid segregation; and regulation of cell cycle process. Located in nucleoplasm.

Recommended Dilution: WB,1:500 - 1:2000 IHC-P,1:50 - 1:200

Synonyms: SORORIN; CDCA5

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 153-252 of human CDCA5 (Q96FF9).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.