

## CAB4302

---

### Product Information

<b>Product SKU:</b> CAB4302	<b>Gene ID:</b> 9456	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b> 42kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 40kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** This gene encodes a member of the homer family of dendritic proteins. Members of this family regulate group 1 metabotropic glutamate receptor function.

**Recommended Dilution:** WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200

**Synonyms:** HOMER; SYN47; Ves-1; HOMER1A; HOMER1B; HOMER1C; HOMER1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 125-354 of human HOMER1 (NP\_004263.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.