

## CAB4321

---

### Product Information

<b>Product SKU:</b>	CAB4321	<b>Gene ID:</b>	9570	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

### Additional Information

<b>Observed MW:</b>	24kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	25kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a trafficking membrane protein which transports proteins among the medial- and trans-Golgi compartments. Due to its chromosomal location and trafficking function, this gene may be involved in familial essential hypertension.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** Bos1; EPM6; GS27; MYOS; GOSR2

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human GOSR2 (NP\_004278.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.