

CAB4392

Product Information

Product SKU: CAB4392	Gene ID: 10247	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 14kDa	Conjugate: Unconjugated
Calculated MW: 14kDa	Isotype: IgG

Immunogen Information

Background: Enables 2-iminobutanoate deaminase activity and mRNA binding activity. Involved in mRNA catabolic process; mRNA destabilization; and organonitrogen compound catabolic process. Located in cytoplasm and nucleus.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: PSP; P14.5; UK114; HRSP12; hp14.5

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-137 of human HRSP12 (NP_005827.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.