

## CAB4398

---

### Product Information

<b>Product SKU:</b> CAB4398	<b>Gene ID:</b> 10286	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b> 25kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 26kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** Involved in mRNA splicing, via spliceosome. Located in centrosome and nuclear speck. Part of U2-type catalytic step 2 spliceosome. Colocalizes with DNA replication factor A complex. Implicated in breast cancer.

**Recommended Dilution:** WB,1:500 - 1:2000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

**Synonyms:** DAM1; SPF27; Snt309; BCAS2

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-225 of human BCAS2 (NP\_005863.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.