

CAB4399

Product Information

Product SKU: CAB4399	Gene ID: 10291	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Mouse,Rat

Additional Information

Observed MW: 130kDa	Conjugate: Unconjugated
Calculated MW: 89kDa	Isotype: IgG

Immunogen Information

Background: This gene encodes a subunit of the splicing factor 3a protein complex. The splicing factor 3a heterotrimer is a component of the mature U2 small nuclear ribonucleoprotein particle (snRNP). U2 small nuclear ribonucleoproteins play a critical role in spliceosome assembly and pre-mRNA splicing.

Recommended Dilution: WB, 1:1000 - 1:2000

Synonyms: PRP21; PRPF21; SAP114; SF3A120; SF3A1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 715-785 of human SF3A1 (Q15459).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.