

CAB4417

Product Information

Product SKU:	CAB4417	Gene ID:	10426	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	105kDa	Conjugate:	Unconjugated
Calculated MW:	104kDa	Isotype:	IgG

Immunogen Information

Background: Enables gamma-tubulin binding activity. Predicted to be involved in meiotic cell cycle; microtubule cytoskeleton organization; and mitotic cell cycle. Located in cytoplasm and microtubule cytoskeleton.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: 104p; ALP6; GCP3; Spc98; SPBC98; Spc98p; Grip104; TUBGCP3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human TUBGCP3 (NP_001273207.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.